


## Cleaning up of $\alpha$ -Methylstyrene by flash chromatography

A. Talamona, BÜCHI Labortechnik AG

**Assessment:** The isolation of  $\alpha$ -Methylstyrene from the crude reaction mixture seems easy. However the reaction mixture contains some dimethyl sulfoxide and the oily wetting agent from the sodium hydride. Both of them are not UV-active and therefore in the TLC invisible.

### TLC of the crude reaction mixture:

TLC on silica gel Si60<sub>F254</sub>, developed in n-hexane/ethyl acetate 90:10, detection UV 254 nm




### Sepacore configuration

Cartridge 40 x 150mm, prepacked with silica gel 60, 40 – 63  $\mu$ m  
2 Pump Modules C-605  
Fraction Collector C-660  
Control-Unit C-620 with SepacoreControl software  
UV Photometer C-635

### Separation conditions

Eluent: n-hexane with 0%, 2% und 15% ethyl acetate, step gradient  
Flow rate: 100ml/min  
Sample: 3 g crude mixture, dissolved in toluene  
Injection volume: 4.5 ml (although the mixture is liquid, the sample is not fully miscible with n-hexane and must be dissolved in toluene and injected directly onto the cartridge).

### Separation



### Fractions

2–10: mainly toluene and wetting agent  
11–12:  $\alpha$ -Methylstyrene  
14–18: Acetophenone (reactant)

### TLC of the separated fractions

TLC on silica gel 60<sub>F254</sub>, developed in n-hexane-ethyl acetate 9:1, detection UV 254 nm

### Recovery

Fraction 11 - 12: 0.4 g colourless oily liquid

